



Product Factsheet  
STA

# Hello and welcome to Sta!

Our vertical and prone stander has been lovingly created to enhance the life of the user. It aims to make standing enjoyable and practical for your child, while simple to achieve.

To gain these benefits, it's important that the Sta is correctly fitted and that's why our customer support team and technical product advisors are here to work with you to find that optimal position. We're here to help.

Wishing you all the very best with your Sta.

The Smirthwaite Team



# Features and Benefits

## Adjustability

Sta is an entry level standing frame for 6-12 year olds. The Sta's anterior pelvic and knees supports are height, width and angle adjustable to achieve midline positioning that helps develop and maintain a child's weight bearing ability. The standing frame is practical too. Sta can be angled with the child in place to achieve the right standing position and as the frame is angled, the table angles too so that a child can continue with their activities.

## Design and Comfort

The Sta has a simple design that appeals to older children wanting something fun yet stylish. Its cushioned chest, hip and leg supports provide extra comfort that encourages longer standing times to help build strength and relieve pressure from sitting. Not only does the Sta deliver physical benefits, but there are also psychological benefits too. Standing allows for greater interaction with family and friends that helps a child improve their social skills.

# Product Specifications

Size	1
Code	21003
User age	6-12 years
Thoracic support (top to foot board)	693-1000
Pelvic support (top to foot board)	478-808
Knee support (top to foot board)	223-343
Foot board height (floor to foot board, excluding sandals)	193
Max width (between thoracic supports)	460
Tray size (depth and width)	390 x 446
Tray angle (°)	0-60
Base frame (length, width & height)	1000 x 660 x 920/1160
Tilt in space (°)	0*-32
Max user weight (kg)	60
Product weight (kg)	28
Price (£)	599

all dimensions in mm unless otherwise stated  
\*upright position



**Smirthwaite**

**16 Wentworth Road, Heathfield**

**Newton Abbot, Devon, UK, TQ12 6TL**

**+44 (0) 01626 835552**

**info@smirthwaite.co.uk**

**[www.smirthwaite.co.uk](http://www.smirthwaite.co.uk)**