



Product Catalog



supporting children and their families for life

seating • changing & showering • therapy • desks & tables • toileting & bathing • sensory • standing

welcome

Since 1986 Smirthwaite Ltd has developed and retained its reputation as one of the leading designers & manufacturers of special needs equipment, both in the United Kingdom and internationally. Strengthening this position, we launched Smirthwaite USA LLC in 2012. In that time, we have dedicated ourselves to finding and providing effective solutions that offer significant benefits to babies, children & young adults in their daily lives. Working closely with children, parents, caregivers and therapists, we have developed a broad range of seating, standing, toileting, bathing, changing and showering products.



Easy to clean

All of our products are easy to clean. Our wooden chairs have antimicrobial properties that ensure surfaces are kept hygienic. This helps to reduce any risk of cross infection and is especially useful in multi-user environments such as nurseries and schools.

All our upholstery is carefully chosen and all contain antibacterial agents that help against the spread of MRSA, E coli and C dif.

Innovation in our design

Special needs equipment is often overly complex, difficult to use and expensive to maintain. At Smirthwaite, we believe in well made, easy to use products that stand the test of time. Through close consultation with leading therapists, parents and users, we ensure our designs offer numerous clinical benefits while maintaining a customer-centric and innovative approach.

All products come with a 2 year warranty as standard.

Inquiries or Orders

You can order our products quickly and easily online via our website below. Or you can call our Customer Services team today on **212 541 2455** or e-mail **info@smirthwaite.com**

www.smirthwaite.com



standing

buddy 2



seating

sunni 3
brookfield 4
hardrock 6
heathfield 8
hip spica / corner chair 11
felix 12
juni 14
zoomi / zoomi + 18



desks & stools

desks 20
wheely stool 21



added extras

sandals 22



toileting & bathing

bath corner chair 23
combi 24



changing & showering

hi riser 26
hi riser showering 28
mobi changer 30
mobi heavy duty changer 32
mobi showering trolley 34
shower stretcher 36



buddy

Supportive and affordable standing

Smirthwaite’s new stander, the Buddy, is an entry-level vertical stander. Designed for younger children (please use technical data for sizing) who require low to moderate levels of support while weight bearing. Its design allows it to be used as a hip stander suitable for targeted therapy for the development of specific areas such as trunk and shoulder girdle stability and head control.

Inclusive and accessible

The Buddy stander is slimline so it can be used in schools, enabling your child to perform numerous activities without feeling isolated from friends. Its color and friendly character also brightens up any room. With the economic cost we are delivering an entry-level stander that is an accessible option to everyone.

Adjustability

The Buddy stander offers numerous height and width adjustments in the chest and pelvic bands as well as additional angle adjustments in the knee supports and sandals. The Buddy size 2 has an additional height adjustment for taller or faster-growing children.

Buddy stander technical data

all dimensions in inches unless otherwise stated

Size	Size 1 Buddy	Size 2 Buddy 2
Item #	21001	21002
Thoracic support (top) to foot board height	18.9 - 24.4	27.6 - 37.0
Pelvic support (top) to foot board height	13.8 - 18.1	19.9 - 29.3
Knee supports (top) to foot board height	7.1 - 13.4	6.3 - 22.0
Foot board to floor height (excluding sandals)	4.7	4.7
Max width (between thoracic supports)	17.7	17.7
Tray size (depth and width)	16.3 x 20.5	16.3 x 20.5
Base frame (length, width and height)	25.6 x 19.6 x 30.9	31.5 x 24.6 x 39.1
Base frame (length, width and height) (high setting)	N/A	31.5 x 24.6 x 43.8
User weight limit (lbs)	88.2	132.3
Product weight (lbs)	30.9	46.3

Buddy and Buddy 2 come complete with tray, chest band, pelvic band, knee supports and sandals. Straps and tray only available in lime.



sunni

School cafeteria chair

Sunni is all about inclusion at meal times at school. It gives children the opportunity to eat and socialize with their friends during school break times instead of eating separately in their own chairs.

Working with Cerebra

It has been designed along side the Charity Cerebra. Families with a child who has a brain condition face challenges every day. Just to learn, play, make friends, enjoy and experience the world can feel difficult, sometimes impossible. But we believe challenges can be over come.

Cerebra's vision is that every family that includes a child with a brain condition will have the chance to discover a better life together.

Their mission is to listen to families that include children with brain conditions. They use what they tell them to inspire the best research and innovation. Cerebra helps them put the knowledge into practice so they can discover a better life together.

Benefits of Sunni

- Designed to fit most cafeteria chairs
- Armrests fold away
- Three point harness
- Wipe clean seat pads
- High backrest for taller children
- Easy to transport and manœuvre
- Can be used at home and in school
- Available in four vibrant colors



Pictured:
Sunni with footrest

Sunni technical data

all dimensions in inches unless otherwise stated

Size	1	1 HB	2	2HB
Seat depth	8.2	8.2	10.1	10.1
Seat width	8.8	8.8	10.8	10.8
Back height	9.9	13.8	11.1	16
Armrest height	5.9	5.9	6.8	6.8
Back angle [°]	3.7	3.7	3.7	3.7
Max user weight (lbs)	33	33	55	55
Product weight (lbs)	6.6	7.7	8.8	9.9

Optional footrest

An optional footrest is available to give support to the legs when sitting at the table.

Colors



brookfield



The Brookfield chair comes complete with moulded back, contoured seat, 4 point pelvic harness, contoured chest harness, back extension and height adjustable armrest as standard.

Increased stability

The Brookfield is a robust chair for children from an early age to young adults who require moderate levels of postural support when seated. The Brookfield has a unique moulded back as standard to promote central alignment while your child is seated.

Greater adjustability

The Brookfield offers seat height and depth adjustability. As an additional feature, the back can also be adjusted in height. This allows even greater flexibility and enables the chair to grow with your child.

Perfect for customization

Each child is unique. The Brookfield has been carefully designed to provide your child with a high level of customization; this means that the needs of your child are brought to the forefront. By using the chair and the accessories together we can offer a highly tailored solution.

Brookfield technical data

all dimensions in inches unless otherwise stated

Item #	Size*	Seat to floor height**	Width	Seat depth	Back height	Max user weight (lbs)	Product weight (lbs)
6571	1LB	10.0 - 13.0	12.8	8.7 - 12.6	14.4 - 19.3	121.3	26
6573	2LB	11.0 - 15.0	14.4	10.8 - 15.2	16.1 - 20.9	143.3	39
6575	3LB	13.6 - 17.5	17.1	14.0 - 18.3	17.7 - 23.2	165.3	51
6577	4LB	16.7 - 22.6	17.1	16.7 - 22.4	19.9 - 25.4	187.4	62

* Each size comes with a back extension as part of the package.

**Height of chair with contoured seat.

brookfield accessories



Pommel

ITEM # PL010 & PL012

Encourages hip abduction



4 mobile legs with 3in casters

ITEM # LG001 - LG004

Casters can be locked



Reversible footboard

ITEM # FB002 - FB004

Provides foot support & prevents unwanted chair movement



Adjustable footrest

ITEM # FR002 - FR004

The sides help to keep feet contained. Easy height adjustability



Folding footrest¹

ITEM # FR015 - FR017

Facilitates standing transfer
Foot plates fold away to the side



Footplate & straps²

ITEM # FP003 - FP005

Ensures that feet are firmly in position



Adjustable height push handle

ITEM # PH018, PH020 & PH022

For optimum push position



Mobile dolly base

ITEM # DL002 - DL004

Easy to manoeuvre with swivel multi-positional casters



Pelvic cushions single or double thickness (pair)³

Provides hip and thigh support
Available in various sizes



Tray handrail half width

ITEM # HR004

Helps to stabilize upper limbs
Additionally can be used to attach toys



Knee blocks

ITEM # KN001

Provides knee position
Encourages hip abduction



Tray

ITEM # TR003 - TR005

Provides trunk support
Additionally can be used for activities.

¹ Must be factory fitted at the time this product is ordered

² Only available in conjunction with a Mobile dolly base, Reversible footboard or Skis with half & full footrest

³ Available as double cushions (pair) size 1+1in, 1.5+1in & 2+1in and singular cushions (pair) size 1in, 1.5in & 2in - please ask for ITEM #'s

hardrock



Robust and durable

The Hardrock is an extremely durable and robust chair that has been specifically designed for children who are prone to rock repetitively. The Hardrock provides moderate levels of postural support and is suitable for both school and home environments.

Tailoring to the child’s needs

It can be a challenge to engage children with complex learning difficulties in educational and everyday activities. The Hardrock tackles this issue and offers a positive solution. It provides a safe environment and one that is conducive to learning.

Posture matters

The Hardrock uses an anti-thrust system and anterior pad to provide secure postural management within clearly defined boundaries. Thoracic supports and stabilizers are available for children who exhibit sideways motion.

Hardrock technical data

all dimensions in inches unless otherwise stated

Item #	Size	Seat to floor height	Width	Seat depth	Back height	Max user weight (lbs)	Product weight (lbs)
8202	2	10.2 - 16.1	13.8	11.8 - 16.3	19.3	132.3	80
8203	3	11.4 - 17.3	16.9	14.6 - 19.1	22.4	154.3	90
8203LB (low back)	3	11.4 - 17.3	16.9	14.6 - 19.1	18.5	154.3	88
8204	4	14.8 - 20.6	18.5	15.7 - 21.7	26.6	198.4	113

Coded for reimbursement A9900 or A9999

hardrock accessories



Tray handrail full width

ITEM # GB010 - GB012

Helps child to stabilize upper limb
Additionally can be used to attach toys



Trunk supports with straps

ITEM # SP010 - SP011

Provides additional trunk support



Pommel

ITEM # PL020 - PL021

Encourages hip abduction



Tray retaining strap set

ITEM # ST056

Secures the tray and keeps the tray
in the desired position



Moulded pelvic support

ITEM # SP017 - SP019

Provides snug contoured support

Hardrock chair comes complete with mobile base, 4 point pelvic belt and either a high or low back (8203LB) slotted upholstery as standard.



Samba style headrest

ITEM # HT014 - HT015

Provides extra support for the head
Can be shaped by adjusting the sides
Available in two sizes:

Size 1 - 12.7 (w) x 4.7 (h)

Size 2 - 15.2 (w) x 4.7 (h)



Footplate & heel huggies

ITEM # FP017 - FP019

Provides extra foot support



heathfield

The perfect start

Heathfield is the perfect starter chair for very young children with sizes up to young adults who require minimal levels of postural support when seated. Its unobtrusive design means that it fits well into nursery, school and home environments.

A natural choice

We believe in wood as a natural material for our chairs. Wood is extremely durable and it holds a contemporary beauty that makes it aesthetically pleasing while being environmentally friendly. Our timber comes from sustainable sources and offers thermal properties that make it warm to the touch.

Meeting your child's needs

The Heathfield and its many accessories work in partnership together to provide support and comfort for your child. We recognize that every child is different and the numerous configurations of the Heathfield will ensure that the chair is tailored to your child's needs.



The Heathfield chair comes complete with a buckle pelvic strap, adjustable armrests and a contoured seat.


Heathfield technical data

all dimensions in inches unless otherwise stated

Item #	Size	Seat to floor* height	Width	Seat depth	Back height	Max user weight (lbs)	Product weight (lbs)
6901	1	7.1 - 9.1	9.4	7.1 - 9.4	9.8	55.1	11
6902	2	7.9 - 11.8	11.0	8.7 - 12.2	11.0	99.2	13
6903	3	9.4 - 14.6	12.8	10.2 - 13.4	12.6	110.2	18
6904	4	11.4 - 15.4	14.4	12.6 - 15.7	14.4	132.3	35
6905	5	12.8 - 16.7	17.1	14.6 - 18.1	16.5	154.3	44

* Height of chair with contoured seat. The height with a flat seat is 1 inch less

heathfield accessories



Skis

ITEM # SK001 - SK005


Provides greater stability to avoid tippy-ness



Ski & full footrest

ITEM # SK023 - SK027

Provides foot support and offers greater stability



Mobile dolly base

ITEM # DL006 - DL010

Easy to move with swivel multi-directional casters



Reversible footboard

ITEM # FB006 - FB010

Foot support that prevents unwanted chair movement



Adjustable footrest

ITEM # FR005 - FR009

Easy to adjust footrest
The sides of the footrest help to keep feet within the chair



Minimum raise mobile base

ITEM # MB005 - MB007

Facilitates standing transfer
Easy to maneuver
Only available with sizes 3 - 5 chairs



4 Point chest harness

ITEM # HN001 - HN004



Pelvic harness

ITEM # HN039 - HN041




heathfield accessories



Pelvic cushions – double thickness (pair)¹

ITEM # CS050 - CS053, CS230-
CS234, CS077 - CS080

Decreases seat width
Available in various sizes



Pelvic cushions – single thickness (pair)²

ITEM # CS005 - CS008, CS014 -
CS017, CS023 - CS026

Provides hip and thigh support
Available in various sizes



Pommel⁵

ITEM # PL001 - PL002, PL013 -
PL014

Encourages hip abduction



Footplate & straps³

ITEM # FP001 - FP005

Ensures feet are firmly in position



Knee blocks^{3,5}

ITEM # KN001

Maintains knee position
Encourages hip abduction
Only available with sizes 3 - 5 chairs



Adjustable height push handle

ITEM # PH010 - PH014

Adjusts to achieve optimum push position. To be used in conjunction with a Mobile dolly base or Minimum raise mobile base



Tray⁴

ITEM # TR001 - TR005

Promotes trunk support
Additionally can be used for activities, playing and feeding



Quick release tray⁶

ITEM # TR101 - TR104

The tray comes in white and only available on chairs sizes 1-4

¹ Available in sizes 1+1in, 1.5+1in & 2+1in

² Available in size 1in, 1.5in & 2in

³ Only available in conjunction with a Mobile dolly base, Reversible footboard or Skis with half & full footrest

⁴ Heathfield chair sizes 1 & 2 must have a suitable stabilizer (footrest, mobile dolly base or skis) when attaching trays

⁵ The Pommel and Knee Blocks cannot be used together

⁶ Heathfield chair sizes 1 & 2 must have a suitable stabiliser e.g. footrest/dolly/skis in order to accommodate tray

corner chair



Maximum adjustability

The Corner Chair offers high level of adjustability in the back, side angle, side pad, headrest and pommel. This means that the chair can be positioned to suit your child's individual needs while also enabling the chair to easily grow with your child. Maximum user weight of 187lbs.

Floor level support

Our Corner Chairs are ideal for young children who need to develop sitting balance. During floor sitting some children are inclined to fall back or sideways if they overstretch or lose concentration. Our Corner Chairs are designed to provide support at the back and sides of children to improve their sitting balance. This encourages social interaction and allows them to engage in floor-based activities without the risk of falling over. The pommel at the front is designed to keep legs abducted to promote a stable posture.

Corner Chair technical data

all dimensions in inches

Item #	Size	Height range	Width	Seat depth
5501	1	16.1-20.4	15.7	11.8
5502	2	20.4-25.9	15.7	11.8
5503	3	21.7 - 27.2	17.7	13.7

hip spica

Support and comfort

Our range of Hip Spica chairs has been specifically designed for young children who are in plaster casts following hip surgery and require supportive seating.

Portable Hip Spica chair

The new design offers high levels of adjustability and creates a lightweight solution. The Portable Hip Spica features a new slimline design and can easily be dismantled so it's simple to transport. It is supplied with a 5-point harness and quick release tray as standard.



Portable Hip Spica chair technical data

all dimensions in inches

Item #	Size	Seat height	Footprint	Tray size	Max user weight lbs	Product weight lbs
5596	1	8.3 - 14.2	22.0 x 26.0	21.7 x 17.3	55	22

felix

simple stackable school chair



Comfortable learning

Standard classroom chairs can be uncomfortable for a whole day of learning but Felix comes with a seat and back cushion all in one, this gives that little extra padding on the seat and back. The Felix chair provides postural support when needed, when a child is comfortable learning becomes more effective.

Simple features

Felix is stackable, an important feature for any classroom at the end of a working day. The cushions are easily removable and wipeable. They come in six vibrant colors and the chair comes in five sizes. Felix also comes flat packed with easy to follow assembly instructions reducing packaging and making it easier for delivery and handling.

Classroom chair

Felix has been designed as an active classroom arm chair. The lateral side arm rests not only provide support while seated but it also facilitates the sit to standing movement and vice versa. The shape of the seat and back promotes good posture and the design of the chair legs provides stability and prevents tipping.

felix

Felix technical data

all dimensions in inches unless otherwise stated

Size	Seat to floor height	Seat Width	Seat depth	Max user weight (kg)	Product weight (kg)
1	10	11 1/2	8 1/2	165	8
2	12	13	10 1/2	176	10.5
3	14	14 1/2	12 1/2	198	14
4	16	16	14 1/2	242	18
5	17 1/2	16	14 1/2	242	22

Colours	Size 1 Item #	Size 2 Item #	Size 3 Item #	Size 4 Item #	Size 5 Item #
pink	8166	8167	8168	8169	8170
teal	8171	8172	8173	8174	8175
lime	8176	8177	8178	8179	8180
cobalt	8181	8182	8183	8184	8185
red	8101	8102	8103	8104	8105
blue	8161	8162	8163	8164	8165



juni



Whether in the classroom, nursery, or at home, good posture when seated creates a positive experience for your child; fine motor skills develop, confidence increases and the ability to engage and learn becomes easier.

The Juni chair has a contemporary appeal that helps to promote social inclusion.

Posture at its foundation

Juni's seat and backrest both provide unique contouring. This promotes better posture, gives higher levels of comfort and helps your child feel more secure and stable in the Juni. It is our vibrant chair for children who require minimal levels of postural support when seated.

Highly adjustable

We understand that your child's needs change during each stage of growth and development. Seat height and depth as well as armrest height can all be easily adjusted. The backrest can be reclined by 5° increments (from 95° to 105°), while it also features an unique built-in lumbar support.

The Juni chair comes with contoured seat and backrest, built-in lumbar support, pelvic strap, height and depth adjustable seat, height adjustable armrests as well as angle adjustable backrest as standard.

Juni technical data

all dimensions in inches unless otherwise stated

Size	1	2	3	4
Seat to floor height	6.3 - 8.7	7.9 - 11.0	9.4 - 13.4	11.8 - 15.7
Width	9.4	11.0	12.8	14.4
Depth	7.1 - 9.4	8.7 - 11.8	10.2 - 14.2	12.6 - 16.5
Back height	10.2	11.4	12.8	15.7
Back angle (°)	95/100/105	95/100/105	95/100/105	95/100/105
Max user weight (lbs)	66.1	88.2	110.2	132.3
Product weight (lbs)	9	13	17	22

Color	Item #'s			
Blue	16031	16032	16033	16034
Lime	16021	16022	16023	16024



juni accessories

We've worked in close collaboration with therapists, parents and children to ensure that our products meet and surpass your expectations. Through your valuable feedback the Juni also has companions: highly adjustable Juni Table, Juni Hi & Lo Boosters and Juni Hi Lo Base to overcome height issues at home, nursery & school.

Juni Table

The Juni table offers two key benefits; the legs are adjustable so your child can work or play at the correct height and the table folds, which allows for easy storage. Perfect companion to your favorite Juni chair.



Juni table technical data

all dimensions in inches unless otherwise stated

Size	1	2
Height	12.2 - 17.7	15.7 - 23.6
Width/length	27.6	30.7
Depth	20.9	23.2

Color	Item #	
Blue top	2431	2432
Lime top	2421	2422

Juni Boosters

Our Hi Booster is perfect for the breakfast bar or for feeding your child without crouching or kneeling. Our Lo Booster allows for sitting your child at the dinner table or giving your child extra height.

Hi & Lo Boosters come complete with a black 5-point harness and are also available with castors*. The Juni Boosters accommodate Juni chair size 1 and 2 only. Please visit www.smirthwaite.com for more details.



Juni Boosters technical data

all dimensions in inches unless otherwise stated

Booster Type	Lo Booster		Hi Booster	
Size of chair on booster	1	2	1	2
Seat to floor height	26.4 - 30.3		34.3 - 38.2	
Footprint open	20.5 x 20.5		20.5 x 25.2	
Folded size	20.5 x 13.4 x 27.6		20.5 x 16.5 x 36.6	

Color	Item #			
Natural	JB411	JB421	JB412	JB422

* Castors with brakes are sold as an accessory (Item # JB001); this will increase the height of the Juni Booster by 2.56 inches

juni accessories



Juni Hi Lo Rolling Base

We understand the challenges of the classroom and home environment. Height adjustment is now easy with our Juni Hi Lo Base. The slimline footprint of the base makes this the ideal solution to fit under most classroom desks. So now, you can use your Juni chair together with the Juni Hi Lo Base by simply attaching it to the plate on top of the base; or use the chair independently where height adjustment's are not needed.

Dolly Base

ITEM # PH101 - PH104
Facilitates manoeuvrability of your Juni chair
Available in 4 sizes



Juni Hi Lo Base technical data

all dimensions in inches unless otherwise stated

Size	1		2	
Dimension	Min	Max	Min	Max
Seat to floor height	6.3	22.4	7.9	22.4
Lower leg length to footrest	6.3	8.7	7.9	11.0
Push handle height	19.7	39.4	19.7	39.4
Base length and width	25.5 x 20.5		25.6 x 20.5	
Max user weight (lbs)	66.1		88.2	
Product weight (lbs)	23.1		23.1	
Color	Item #'s			
Natural Footrest	HL141		HL142	

juni accessories



Skis

ITEM # SK101 - SK104

Helps to stabilize and prevent chair movements



Mobile skis & footrest

ITEM # MK101 - MK104

Gives clearance to facilitate standing transfer. Includes a sliding footrest & push handle.



Ski & footboard¹

ITEM # SK111 - SK114

Provides foot support and offers greater stability



Footrest

ITEM # FR051 - FR054

Height and angle adjustable footrest. Can be attached to the Juni or used separately.



4-point harness

ITEM # HN231 - HN234

Provides additional trunk support



Quick release tray

ITEM # TR101 - TR104

White quick release tray



Side support cushions

Provides hip and thigh support
Available in 1 or 2 inch sizes. Sold as a pair.

ITEM # - 1 Inch CS531 & CS532

ITEM # - 2 Inch CS631 & CS632



Pommel

ITEM # PL131 & PL132

Encourages hip abduction. Can be removed with a hand tool.



zoomi

Benefits of Zoomi

- Economical seating solution
- Grows with your child
- Highly adjustable
- Tailored to the needs of your child
- Contemporary design
- Easy to clean
- Stable base with slimline footprint
- Can be used at home and in school

Our standard Zoomi

The standard Zoomi is ideally suited to the classroom and home. Both the seat and footrest can be easily adjusted to the appropriate height and depth, allowing the chair to accommodate your child's needs as they grow.



Standard features

- Frame
- Low backrest
- Flat laminated seat
- Laminated footboard
- Pelvic strap

Zoomi technical data

all dimensions in inches unless otherwise stated

Size		1
Seat to floor height range		4.3 - 24.2
Seat width		17.7
Seat depth		21.7
Back height		10.2
Back angle (°)		95/100/105
Frame height		31.9
Max user weight (lbs)		187.4
Product weight (lbs)		35
Colour	Zoomi standard	Zoomi +
Blue	13031	13032

Perfectly customizable

With multiple options to cater for a variety of users, Zoomi can be uniquely configured to meet the needs of your child. We have put together two simple options:

- Standard Zoomi
- Zoomi+

zoomi+

One size fits all

Zoomi's innovative and unique seating system has been designed to grow with your child, covering an extensive age range, all the way up to young adults.

Zoomi+

Following collaboration with a number of therapists, we have developed Zoomi+. Complete with low back padded cover, armrest, quick release tray, 5-point harness, flat padded seat, slotted footrest and mobile skis, it offers the ideal solution for children who require a little extra postural support.

Standard features

- Frame
- Low back padded cover
- Armrest
- Quick release tray
- 5-point harness
- Flat padded seat
- Slotted footboard
- Mobile skis




Armrest

ITEM # AR019

Provides arm support

Used for the Quick release tray



Quick release tray

ITEM # TR104

Helps to stabilize trunk

Available in white only

Must be used with armrest

desks

study desk



activity table



Perfect for school and home

Our range of desks have been designed to fit in to the school and home environment. The contemporary beauty of wood and the quality of the design enables our desks to stand the test of time.

Study Desk

The Study Desk has been cleverly designed to allow your child to take part in various activities. The desk also contains ample storage space, so books and toys can be hidden away after they have been used. The Study Desk can be used with a variety of chairs, including the Heathfield, Brookfield and Juni.

Study Desk technical data

all dimensions in inches

Item #	Size	Height range	Width	Depth	Max user weight (lbs)	Product weight (lbs)
6967	1	13.6 - 20.5	17.7	23.6	176	42
6968	2	20.5 - 32.3	21.7	29.5	176	42

Practical and easy to store

Our Activity Table is an ideal height-adjustable table for younger children at nursery, school and home. The steel legs can be adjusted in height to ensure good seating posture of your child and can easily be folded for storage.

Activity Table technical data

all dimensions in inches

Item #	Size	Height range	Length	Width	Max user weight (lbs)	Product weight (lbs)
2251	1	15.2 - 23.4	29.5	23.6	176	26
2252	2	21.9 - 30.1	29.5	23.6	176	26



wheely stools



Therapist’s favorite

The Adjustable Height Stool is our most popular wheely stool. It has been designed specifically for working at lower levels with your child, student or client. The stool is designed for sitting along the narrower side. Available in green and blue with a height range of 13.4 - 18.1in. One size only.

Adjustable Height Wheely Stool technical data

all dimensions in inches unless otherwise stated

Item #	Size	Color	Width	Height range	Max user weight (lbs)	Product weight (lbs)
4182	Standard	Green & Blue	17.7	13.4 - 18.1	220.5	16
4193	Heavy Duty	Green & Blue	17.7	13.4 - 18.1	330.6	17.6



sandals

Feet support

Children's feet also require additional support to help maintain the position of their hips. Feet support enhances good balance and enables freedom of movement in the upper body.



Ferndale sandals (pair)

ITEM # SL001 - SL004

Available in 4 sizes



Cocoon sandals (pair) – velcro¹

ITEM # (standard) SL005,
SL007, SL009, SL011, SL013



Heel huggies (pair)

ITEM # HH001 - HH005

Available in 5 sizes

Ferndale Sandals Fitting Measurements

all dimensions in inches

Item #	Size	Length	Width
SL001	1	7	3.9
SL002	2	8.2	4
SL003	3	9.4	4.2
SL004	4	11	4.3

Cocoon Sandals Fitting Measurements

all dimensions in inches

Item #	Size	Length	Width (heel to toe)
SL005	0	5.5	1.9-2.9
SL007	1	7.4	2.7-3.9
SL009	2	8.2	3.3-4.4
SL011	3	10.6	3.5-4.5
SL013	4	11.2	3.4-4.9

Heel Huggies Fitting Measurements

all dimensions in inches

Item #	Size	Depth	Max Circum	Min Circum
HH001	1	1.9	6.2	5.1
HH002	2	1.9	7.8	6.6
HH003	3	2.9	9.4	8.2
HH004	4	2.9	11	9.8
HH005	5	3.5	12.2	11

bath corner chair



The Bath Corner Chair comes complete with velcro pelvic strap, adjustable pommel, telescopic back rail and telescopic handrail as standard.

More support

Our Bath Corner Chair is designed for children who require a little extra help with sitting balance when bathing. This is achieved by its angled back and hand rail, both of which provide additional security.

Adjustable

Like the Bath Chair, our Bath Corner Chair provides numerous levels of adjustment. The pommel and pelvic strap can be adjusted to suit each individual child's needs.


Compact

Lightweight and easy to clean, the Bath Corner Chair has a perfectly compact design. This allows it to sit low in the bath enabling your child to fully extend their legs in complete comfort. Its size also enables it to fit a wide range of baths.

Bath Corner Chair technical data

all dimensions in inches unless otherwise stated

Size	1	2
Item #	7511	7512
Seat height	3.0	3.0
Seat width	13.4	13.4
Seat depth	11.0	13.0
Back height	13.8	18.1
4 Pommel settings (distance to back of pommel)	4.3 / 5.1 / 6.3 / 7.5	5.5 / 6.7 / 7.9 / 9.1
Max user weight (lbs)	66.1	77.2
Product weight (lbs)	6.6	8.8



combi



Progressive toileting

Our Combi toileting chair is a great starter chair. Its design enables it to be used as a supportive chair while at the same time providing numerous features to assist your child with the everyday toileting needs.

Providing solutions

The Combi toileting chair’s design also incorporates a potty for younger children. This can easily be removed to provide a commode seat for children who have progressed their toilet training.

Adjustability

Good levels of adjustability are offered with the Combi toileting chair so that a wide range of needs can be met. Seat height and depth can be adjusted to suit the needs of your child.

The Combi Chair comes complete with potty, plastic commode seat, upholstered seat and velcro pelvic strap as standard.

Combi technical data

all dimensions in inches unless otherwise stated

Size	1	2	3	4
Item #	7551	7552	7553	7554
Seat to floor height*	8.5 - 11.4	10.8 - 14.8	11.8 - 15.7	13.4 - 17.7
Seat width	11.4	12.8	14.4	17.1
Seat depth	8.1 - 11.4	9.6 - 14.0	12.2 - 15.2	13.8 - 17.3
Back height	12.8	14.6	15.7	19.7
Max user weight (lbs)	77.2	99.2	121.3	143.3
Product weight (lbs)	24.3	28.7	30.9	33.1

*Seat height to floor is to the top of upholstered seat (deduct 1in for the height of the Potty and deduct 1.6in for the height of the Commode seat)

combi accessories



Skis & full footrest¹

ITEM # SK037 - SK040

Provides foot support and offers greater stability



Reversible footboard¹

ITEM # FB023 - FB026

Provides foot support
Prevents unwanted chair movement



Independent adjustable footrest¹

ITEM # 4187

Portable, height and angle adjustable
Solid support for the feet



Skis

ITEM # SK051 - SK054

Provides greater stability



Footrest & straps¹

ITEM # FP002 - FP005

Ensures feet are firmly in position



Splash guard

ITEM # SG001 - SG002

Assists in leg abduction and keeps toileting contained
Requires Toggle handrail



Pelvic cushions single or double thickness (pair)²

Provides hip and thigh support
Available in various sizes



Toggle handrail

ITEM # HR014 - HR017

Allows child to hold and maintain balance and to sit with minimal support



Tray

ITEM # TR003 - TR006

Provides trunk support
Additionally can be used for activities such as playing

¹ Only available in conjunction with a Mobile dolly base, Reversible footboard or Skis with half & full footrest

² Available as double cushions (pair) size 1+1in, 1.5+1in & 2+1in and singular cushions (pair) size 1in, 1.5in & 2in - please ask for ITEM #'s

hi-riser

(powder coated for dry rooms and stainless steel for bath rooms)



Robust yet slimline

Hi-Riser is an electrically operated wall-mounted changing table that is suitable for children and adults of up to 441lbs. Despite its robust structure, the Hi-Riser has a slimline design that ensures it fits easily into most environments.

Easy transfers


The Hi-Riser has been designed so that it can be used with most mobile or ceiling hoists. The unique safety side rail stows away, streamlining access for transfers from a wheelchair, shower and toileting chair. The Hi-Riser also benefits from an impressive travel range of 24in in height.

Providing options

We understand that every need is different, which is why the Hi-riser comes in seven different sizes that cater from young children through to adults. The Hi-Riser comes in a cream powder coated frame for dry-rooms or stainless steel frame for bath-rooms (not suitable for direct water contact).


hi-riser

key features



- Electrically operated
- Wall and floor mounted
- Cream powder coated frame for dry rooms
- Stainless steel for bath rooms
- Weight capacity of 441lbs
- Impressive travel range of 24"
- Folds away to save on space
- Easy to clean padded surfaces
- Wall mounting kit included
- Stowaway side rail

Mild Steel Powder Coated	Stainless Steel	Size	User Weight	Table Dimensions			Travel Range	Wall Frame Width	Wall Frame Depth	Wall Frame Height	Safety Side Height	Stored Width	Shipping Weight
Item #	Item #		lbs	L in	W in	M in	H in	A in	B in	C in	Y in	S in	lbs
3043	3053	3	441	55.1	23.6	36.4	15.7 - 39.4	27.6	4.3	61.0	6.5	12	185.2
3047	3057	7	441	70.9	25.6	38.4	15.7 - 39.4	27.6	4.3	61.0	6.5	12	202.8



hi-riser shower stretcher



Two in one

The Hi-Riser shower stretcher is an easy to use two-in-one wall and floor mounted showering and changing table. Its robust yet slimline design makes it suitable for children and adults of up to 330lbs. Its versatility means it offers the ideal solution for swimming, hydrotherapy and bathroom environments.

Comfortable & hygienic

The Hi-Riser Shower Stretcher has been designed with reinforced PVC, making it highly supportive yet soft and comfortable to lie on. Its stainless steel frame facilitates easy cleaning, making it a perfect choice for multi-user environments.

Quick drainage

The Hi-Riser Shower Stretcher is supplied with flexible pipe work to use with a suitable floor drain or, alternatively, it can be plumbed in. Water can easily be drained away using the drainage pipe, which can also be detached for cleaning purposes.


hi-riser shower stretcher

key features



- Electrically operated – handset works safely in wet rooms
- Changing and showering in one
- Drainage pipe & fasteners included
- Wall and floor mounted
- Stainless steel construction
- Weight capacity of 331lbs
- Travel range of 22.4" in height
- Easy to clean padded surfaces
- Unique stowaway side rail
- Folds away to save on space
- Only protrudes 11.8" from the wall when stored

Stainless Steel	Size	User Weight	Table Dimensions			Travel Range	Wall Frame Width	Wall Frame Depth	Wall Frame Height	Safety Side Height	Stored Width	Shipping Weight
Item #		lbs	L in	W in	M in	H in	A in	B in	C in	Y in	S in	lbs
3067	2	331	70.9	25.6	43.7	15.4 - 37.8	27.6	4.3	61	6.5	11.8	181



mobi changer



This image shows the product in a stainless steel option

Room to room

The Mobi Changer is a fully mobile changing table that requires no installation or attachments. Its compact and slimline design means it can be easily maneuvered from room to room, offering an extremely cost effective option for multi-transfer or multi-user environments.

Height adjustable

The electrically operated Mobi Changer provides an excellent travel range of 37.4in in height and can be lowered and elevated at the touch of a button. Safe and suitable working heights can be easily achieved which helps prevent the risk of back injury to the caregiver.

Easy transfers

The Mobi Changer has been designed to ensure sufficient clearance for a wide number of hoists. Its mobility and unique stowaway side rails enable safe and easy handling. Whether you are transferring from a bed, wheelchair or toileting chair, the Mobi Changer is the perfect solution.

mobi changer

key features




This image shows the product in a powder coated option



- Mobile and easy to maneuver
- No installation required
- Cream powder coated frame for dry rooms
- Stainless steel for wet rooms
- Available in 2 sizes
- Length of table available 55.1in or 70.9in
- Weight capacity of 500.4lbs
- Travel range of 37.4in in height
- Easy to clean padded surfaces
- Optional angled back
- Unique stowaway side rails

Mild Steel Powder Coated	Stainless Steel	Size	User Weight	Table Dimensions		Travel Range	Safety Sides		Safety Sides Operational Space	Castor Centres	Nominal Clearances		Shipping Weights
Item #	Item #		lbs	L in	W in	H in	Width X in	Height Y in	Z in	P in	Q in	G in	lbs
3121	3141	1	500.4	55.1	29.5	18.5 - 37.4	33.1	6.5	7.9	42.9	27.6	7.5	189.6
3123	3143	3	500.4	70.9	29.5	18.5 - 37.4	33.1	6.5	7.9	42.9	27.6	7.5	198.4



mobi heavy duty changer



Meeting bariatric needs

The Mobi Heavy Duty Changer has been designed to meet the needs of bariatric users. Its powerful 12K Newton motor, coupled with the extra wide mattress platform, makes it the ideal solution.

Maneuverability

The benefit of the Mobi range of products is that they do not require any wall mounting or floor attachments. They are easy to maneuver which means that the user can be transported from room to room with no trouble (door opening has to be 43.3in or more).

Height adjustable

The Mobi Heavy Duty Changer is electronically operated, which makes adjusting the height of the table simple, this is beneficial for both the caregiver and user.


mobi heavy duty changer

key features



- Mobile and easy to maneuver
- No installation required
- Electronically operated
- Cream powder coated frame for dry rooms only
- Length of table available in 78.7in
- Weight capacity of 1000lbs
- Travel range of 36.2in in height
- Easy to clean padded surfaces
- Unique stowaway side rails
- Wall mounted charger and 2 batteries come as standard

Mild Steel Powder Coated	Size	User Weight	Table Dimensions		Travel Range	Safety Sides		Safety Sides Operational Space	Castor Centres		Nominal Clearances	Shipping Weights
Item #		lbs	L in	W in	H in	Width X in	Height Y in	Z in	P in	Q in	G in	lbs
3161	1	1000	78.7	37.8	21.7 - 36.2	41.3	6.5	7.9	59.1	31.5	5.9	220.5



mobi shower trolley



Safe showering

For many children and adults, conventional methods of bathing and showering can be difficult and even hazardous. The Mobi Shower Trolley is a powered, height-adjustable shower trolley specifically designed to enable the user to be safely and easily washed or showered by their caregiver.

Reduces back injury


The Mobi Shower Trolley offers a range of height adjustments to enable the caregiver to shower and change the user at a safe and suitable working height. This reduces the risk of back injury for the caregiver.

Stowaway side rails

The Mobi Shower Trolley has been designed with two unique safety side rails that can be easily stowed away beneath the trolley. This allows maximum clearance for safe, easy transfer and use.


mobi shower trolley

key features



- Mobile and easy to maneuver
- No installation required
- Stainless steel frame
- Drainage channel and pipe
- A soft comfortable mattress and wedge shaped pillow
- Available in 1 size
- Length of table is 70.9in
- Weight capacity of 485lbs
- Travel range of 36.4in in height
- Easy to clean padded surfaces
- Unique stowaway side rails

Stainless Steel	Size	User Weight	Table Dimensions		Travel Range	Safety Sides		Safety Sides Operational Space	Caster Centres		Nominal Clearances	Shipping Weights
Item #		lbs	L in	W in	H in	Width X in	Height Y in	Z in	P in	Q in	G in	lbs
3151	1	485.0	70.9	29.5	18.3 - 36.4	33.1	7.7	7.9	42.9	27.6	7.5	189.6



shower stretcher fixed height



Economical solution

The Shower Stretcher Fixed Height is our entry-level product and an economical washing and showering solution for the user. This durable and slimline product represents real value for money.

Makes bathing easier

The Shower Stretcher Fixed Height is designed to allow the washing or showering of the user while in a horizontal position. It can be fixed directly above a bath making it the perfect solution for the home.

Added headrest

The Shower Stretcher Fixed Height, along with all of our shower stretchers, can be enhanced with the optional integral vinyl headrest. The headrest features a paddle lift adjustment, which allows it to be locked securely into the desired position.

shower stretcher fixed height


key features



- No floor attachments required
- Ideal for a wet room environment
- Stainless steel frame with powder coated finish
- Permanently fixed at the desired height
- Available in 2 sizes
- Length of stretcher available 55.1in & 74.8in
- Weight capacity of 440.9lbs
- Folds away to save on space
- Wipe clean surfaces

Stainless Steel Powder Coated		Size	User Weight lbs	Table Dimensions			Travel Range	Wall Frame Width A in	Wall Frame Depth B in	Stored Width S in	Product Weight lbs
Without Headrest Item #	With Headrest Item #			L in	W in	M in					

3171	3271	1	440.9	55.1	25.6	39.4	Fixed	26.2	3.1	9.4in	165.3
3176	3276	6	440.9	74.8	25.6	39.4	Fixed	26.2	3.1	9.4in	176.4



seating standing toileting & bathing changing therapy showering sensory

a helping hand around the World



telephone
212 541 2455

fax
917 368 8005

email
info@smirthwaite.com

Smirthwaite USA LLC
405 Lexington Avenue, 26th Floor
New York NY10174

www.smirthwaite.com

Smirthwaite reserves the right to modify product specifications at anytime without prior notification.
All details shown are correct at time of print. Smirthwaite® is a registered trade mark of Smirthwaite Limited.

