

New
Product!



STAR Product Factsheet

Hello and welcome to Star!

Our vertical and prone stander has been lovingly created to enhance the life of the user. It aims to make standing enjoyable and practical for your child, while simple to achieve.

To gain these benefits, it's important that the Star is correctly fitted and that's why our customer support team and technical product advisors are here to work with you to find that optimal position. We're here to help.

Wishing you all the very best with your Star.

The Smirthwaite Team



Features and Benefits

Adjustability

Star is an entry level standing frame for 6-12 year olds. The Star's anterior pelvic and knee supports are height, width and angle adjustable to achieve midline positioning that helps develop and maintain a child's weight bearing ability. The standing frame is practical too. Star can be angled with the child in place to achieve the right standing position and as the frame is angled, the table angles too so that a child can continue with their activities.

Design and Comfort

The Star has a simple design that appeals to older children wanting something fun yet stylish. Its cushioned chest, hip and leg supports provide extra comfort that encourages longer standing times to help build strength and relieve pressure from sitting. Not only does the Star deliver physical benefits, but there are also psychological benefits too. Standing allows for greater interaction with family and friends that helps a child improve their social skills.



Product Specifications

Size	1
Item	21003
Max user weight (lbs)	132
Min/max user height (approx)	43.3 - 59
Thoracic support (top to foot board)	27-39
Pelvic support (top to foot board)	18.8-31.8
Knee support (top to foot board)	8.7-13.5
Foot board height (floor to foot board, excluding sandals)	7.5
Max width (between thoracic supports)	18.1
Tray size (depth and width)	15.3 x 17.5
Tray angle (°)	0-60
Base frame (length, width & height)	39.3 x 25.9 x 36.2/45.66
Tilt in space (°)	0*-32
Product weight (lbs)	61.7

all dimensions in inches unless otherwise stated
*upright position

Smirthwaite USA LLC

405 Lexington Avenue, 26th Floor

New York

NY10174

Tel: 212 541 2455

info@smirthwaite.com

www.smirthwaite.com