



Postural Seating

standing • seating • desks • stools
toileting • bathing • changing • showering

supporting children and
their families **for life**

smirthwaite

Together we
transform lives

welcome

Since 1986, Smirthwaite Ltd has developed and retained its reputation as one of the leading designers and manufacturers of pediatric equipment, both in the United Kingdom and internationally. Strengthening this position, we launched Smirthwaite USA LLC in 2012. In that time, we have dedicated ourselves to finding and providing effective solutions that offer significant benefits to babies, children and young adults in their daily lives. Working closely with children, parents, caregivers and therapists, we have developed a comprehensive and bespoke range of seating, standing, toileting, bathing, changing and showering products.

Our mission

Is to work together to enhance the lives of people with additional needs using lovingly created specialist equipment to help each person shine.

Lovingly created

Respected throughout the industry for our technical expertise, we play an active role in helping to shape standards and best practice. We believe in well-made pediatric products that are easy to use and, importantly, can stand the test of time. For added convenience, we have produced product user guides and videos which are available to download.

Our durable products are available in a variety of amazing colours to brighten up any child's room (or to simply to allow

them to choose and enjoy their favourite colour chair, for example).

Our bespoke products come in a variety of amazing colours to brighten up any child's room (or to simply to allow them to choose and enjoy their favourite colour chair, for example).

We have a very close relationship with therapists, children and parents which helps us to establish exactly what is required from our products.

Innovation in our design

Pediatric equipment is often overly complex, difficult to use and expensive to maintain. At Smirthwaite, we believe in well made, easy to use products that stand the test of time. Through close

consultation with leading therapists, parents and users, we ensure our designs offer numerous clinical benefits while maintaining a customer-centric and innovative approach. All products come with a 2 year warranty as standard.

Easy to clean

All of our products are easy to clean. Our wooden chairs have antimicrobial properties that ensure surfaces are kept hygienic. This helps to reduce any risk of cross infection and is especially useful in multi-user environments such as nurseries and schools. All our upholstery is carefully chosen and all contain antibacterial agents that help against the spread of MRSA, E coli and C dif.

supporting children & their families for life

For home or school

Whether for home or school, good posture is at the heart of our special needs chairs. Good posture is crucial in so many ways and our paediatric seating range offers some huge benefits:

- Correctly fitted special needs chairs help children focus on their school work or tasks at home
- Helps children to join in with their friends
- By easing pressure on the digestive system, there is less risk of pain and discomfort when eating
- By reducing the stress and strain of poor posture, children can physically and mentally relax

All of the these can lead to a happier, more confident child.

Every child is unique and special and many of our products can be tailored to fit specific needs. Don't be afraid to ask us about customising our products for your child.

Inquiries or Orders

You can order our products swiftly and easily online via our website below; or call our friendly Customer Services team on

212 541 2455
info@smirthwaite.com

seating contents

asiento	4
sunni	7
brookfield	8
hardrock	12
heathfield	14
corner chair / hip spica	17
felix	18
juni	19
zoomi / zoomi +	22
colours	24

*DISCLAIMER - Please note there may be an extended lead time on larger orders that would require transportation on a container versus being able to ship via UPS from the UK.



Asiento

NEW!

Developed alongside input from clinical professionals, the Asiento is a strong and highly adjustable chair suitable for children from 1 year old to young adults with moderate to complex needs. Created to improve postural support and symmetry, the Asiento delivers many clinical benefits while being functional too.

It's durable

Accessories are attached to the chair's strong central column making the Asiento robust and sturdy.

Dynamic back

The dynamic back is great for those who have dystonic body movements. The backrest will move with the user and then gently return them to the correct sitting position.

Easy to use tray

The universal tray is easily attached to the armrests. The tray adjusts with the armrests to allow for use on all seat widths and also depths, providing the best fit possible for all sized users.

Quiet electric hi-lo base

With the electric hi-lo base, users can smoothly raise and lower the chair to the desired height with a user-operated remote control. It's quiet too, helping those who are sensitive to sound.


Closer contact

The open design of the chassis and pushbike handles allows carers to get close to the backrest. This is especially useful for upper-limb therapy, crossing midline, administering care and calming the user after an epileptic fit or when anxious.

designed specifically for children, teenagers and young adults with moderate to complex needs

New deluxe colour range

The following colour options apply to the Asiento:



L Pressure care relief technology in partnership with Levitex has been applied to the Asiento.



Standard key features



Bodypoint harness



Adjustable side supports



Dynamic back & lap strap



Flip-up footrest & multi-positional leg rests

Asiento technical data

all dimensions in inches unless otherwise stated

	Size 1	Size 2	Size 3
Code	1011	1012	1013
Max user weight	176	176	265
Seat depth (min-max)	7.4 - 11.8	10.2 - 14.2	14.4 - 23.6
Seat width	11.2	16.1	18.9
Back rest height	13.8 - 17	16.5 - 20.1	19.3 - 23.2
Seat height from floor	11.4 - 22.2	11.4 - 22.2	19.3 - 23.2
Lower leg length (seat to footplate)	8.3 - 12.6	8.3 - 12.6	12.2 - 22
Armrest height (min-max from seat)	4.3 - 12.8	4.3 - 12.8	4.3 - 12.8
Lateral supports travel adjustment	2	2.6	2.6
Hip guide travel adjustment	2.8	2.8	2.8
Adduction pad travel adjustment	1.9	1.9	1.9
Chest width	9.8 - 13.8	9.8 - 15.4	12.2 - 17.3
Armrest angle (forwards/backwards)	90°/90°	90°/90°	90°/90°
Back angle (forwards/backwards)	10°/32°	10°/32°	10°/32°
Footrest angle (forwards/backwards)	90°/20°	90°/20°	90°/20°
Footplate Angle (plantarflexion/dorsiflexion)	45°/20°	45°/25°	45°/25°
Abduction angle	23°	23°	23°
Adduction angle	23°	23°	23°
Tilt in space (forwards/backwards)	10°/20°	10°/20°	10°/20°
Tray size	26 x 17.3	26 x 17.3	26 x 17.3
Seat width between thighs	9 - 14.6	10.2 - 14.6	13.1 - 17.5



enhancing lives...

sunni

School cafeteria chair

Sunni is all about inclusion at meal times at school. It gives children the opportunity to eat and socialize with their friends during school break times instead of eating separately in their own chairs.

Working with Cerebra

It has been designed along side the Charity Cerebra. Families with a child who has a brain condition face challenges every day. Just to learn, play, make friends, enjoy and experience the world can feel difficult, sometimes impossible. But we believe challenges can be over come.

Cerebra's vision is that every family that includes a child with a brain condition will have the chance to discover a better life together.

Their mission is to listen to families that include children with brain conditions. They use what they tell them to inspire the best research and innovation. Cerebra helps them put the knowledge into practice so they can discover a better life together.

Benefits of Sunni

- Designed to fit most cafeteria chairs
- Armrests fold away
- Three point harness
- Wipe clean seat pads
- High backrest for taller children
- Easy to transport and maneuver
- Can be used at home and in school
- Available in four vibrant colors

Sunni technical data

all dimensions in inches unless otherwise stated

Size	1	1HB	2	2HB
Seat depth	8.2	8.2	10.1	10.1
Seat width	8.8	8.8	10.8	10.8
Back height	9.9	13.8	11.1	16
Armrest height	5.9	5.9	6.8	6.8
Back angle (°)	3.7	3.7	3.7	3.7
Max user weight (lbs)	33	33	55	55
Product weight (lbs)	6.6	7.7	8.8	9.9



Pictured:
Sunni with footrest

Optional footrest

An optional footrest is available to give support to the legs when sitting at the table.

Colors



brookfield



We heard your feedback and have made necessary adjustments to meet your needs. Our recent brookfield and hardrock specifications offer better value and are more cost-effective with quicker lead times, without compromising quality. Expect less fuss and fewer accessories, making it easier to use. Take a look...

New brookfield specifications

- Low, moulded back
- Contoured seat with buckle strap
- Includes standard 2 point strap
- Ocean blue

Increased stability

The Brookfield is a robust chair for children from an early age to young adults who require moderate levels of postural support when seated. The Brookfield has a unique moulded back as standard to promote central alignment while your child is seated.

Greater adjustability

The Brookfield offers seat height and depth adjustability. As an additional feature, the back can also be adjusted in height. This allows even greater flexibility and enables the chair to grow with your child.

Perfect for customization

Each child is unique. The Brookfield has been carefully designed to provide your child with a high level of customization; this means that the needs of your child are brought to the forefront. By using the chair and the accessories together we can offer a highly tailored solution.

Brookfield technical data

all dimensions in inches unless otherwise stated

Item #	Size*	Seat to floor	Width	Seat depth	Back height	Max user weight (lbs)
6571	1LB	10.0 - 13.0	12.8	8.7 - 12.6	14.4 - 19.3	121.3
6573	2LB	11.0 - 15.0	14.4	10.8 - 15.2	16.1 - 20.9	143.3
6575	3LB	13.6 - 17.5	17.1	14.0 - 18.3	17.7 - 23.2	165.3
6577	4LB	16.7 - 22.6	17.1	16.7 - 22.4	19.9 - 25.4	187.4

**Height of chair with contoured seat.

brookfield accessories




Pommel
ITEM # PLO10 & PLO12
Encourages hip abduction



4 mobile legs with 3in casters
ITEM # LG001 - LG004
Casters can be locked



Reversible footboard
ITEM # FB002 - FB004
Provides foot support & prevents unwanted chair movement



Adjustable footrest
ITEM # FR002 - FR004
The sides help to keep feet contained. Easy height adjustability



Folding footrest¹
ITEM # FR015 - FR017
Facilitates standing transfer
Foot plates fold away to the side



Footplate & straps²
ITEM # FP003 - FP005
Ensures that feet are firmly in position



Adjustable height push handle
ITEM # PH018, PH020 & PH022
For optimum push position



Mobile dolly base
ITEM # DL002 - DL004
Easy to manoeuvre with swivel multi-positional casters



Pelvic cushions single or double thickness (pair)³
Provides hip and thigh support
Available in various sizes



Headrest or back extension
ITEM # HR004



Knee blocks
ITEM # KN001
Provides knee position
Encourages hip abduction




Tray
ITEM # TR003 - TR005
Provides trunk support
Additionally can be used for activities.

1 Must be factory fitted at the time this product is ordered
2 Only available in conjunction with a Mobile dolly base, Reversible footboard or Skis with half & full footrest
3 Available as double cushions (pair) size 1+1in, 1.5+1in & 2+1in and singular cushions (pair) size 1in, 1.5in & 2in - please ask for ITEM #'


harness measurements

brookfield, heathfield & combi




Teenage Harness
all dimensions in inches

Size	1	3
A	12.2	14.1
B	11	13.3
C	7	9
D	5.9	6.2




Zipped Jacket
all dimensions in inches

Size	1	2	3
A	2.7	3.5	11
B	9.4	11.4	12.9
C	10	12.5	13.9
D	2.7	3.5	4.3
E	1.9	1.9	2.1




Butterfly Harness
all dimensions in inches

Size	One Size
A	2.7
B	9.4



Bib Harness (padded)
all dimensions in inches


Size	1	2	3
A	4.3	5.5	6.8
B	9	10.4	12.5
C	8.6	11	12.2
D	3.5	5.1	6.1
E	2.3	2.5	2.9



Waitcoat Harness
all dimensions in inches

Size	1	2	3
A	1.9	3.9	5.1
B	5.9	7.8	10
C	12.2	15.3	8.6
D	3.9	5.9	8.6
E	1.9	1.9	2.5
F	24.8	30.7	38.5

F = circumference front to back



Groin Harness
all dimensions in inches

Size	1	2	3
A	9.4	10.6	12.5
B	6.8	9.8	11
C	14.7	18.5	23.6
D	2.5	3.1	3.5
Weight	330	440	661




4-Point Harness
all dimensions in inches

Size	1	2	3
A	9	11	11.4
B	11	11	11.8


harness measurements

samba & hardrock




Soft Cross Harness
all dimensions in inches

Size	1	2	3	4
A	18.1	20	22.4	24.8
B	1.9	1.9	1.9	1.9



Groin Harness
all dimensions in inches

Size	1	2	3	4
A	8.6	8.6	8.6	11
B	6.2	6.2	6.2	8.6
C	6.2	6.2	6.2	6.2
D	1.9	1.9	1.9	2.7
E	10.6	10.6	10.6	10.6
F	1.9	1.9	1.9	2.3



Mid-line support harness with sheepskin pad

Size	1	2	3	4
A (cm)	16	16	16	20
B (cm)	12	12	12	12
C (cm)	7-17	7-17	7-17	7-17

Smirthwaite harnesses offer comfort whilst maintaining secure positioning in all our seating systems. A range of colours and fabrics are available and 'made to measure' solutions and modifications can be made upon request.

Designed for use on Smirthwaite products only.



hardrock

hardrock accessories



New hardrock specifications

- Lowback
- Fully slotted
- Standard 2 point lap strap
- Anti thrust cushions
- Ski transporter base
- Blue

Robust and durable

The Hardrock is an extremely durable and robust chair that has been specifically designed for children who are prone to rock repetitively. The Hardrock provides moderate levels of postural support and is suitable for both school and home environments.

Tailoring to the child's needs

It can be a challenge to engage children with complex learning difficulties in educational and everyday activities. The Hardrock tackles this issue and offers a positive solution. It provides a safe environment and one that is conducive to learning.

Posture matters

The Hardrock uses an anti-thrust system and anterior pad to provide secure postural management within clearly defined boundaries. Thoracic supports and stabilizers are available for children who exhibit sideways motion.

Hardrock technical data

all dimensions in inches unless otherwise stated

Item #	Size	Seat to floor height	Width	Seat depth	Back height	Max user weight (lbs)	Product weight (lbs)
8202	2	10.2-16.1	13.8	11.8-16.3	19.3	132.3	80
8203	3	11.4-17.3	16.9	14.6-19.1	22.4	154.3	90
8203LB (low back)	3	11.4-17.3	16.9	14.6-19.1	18.5	154.3	88
8204	4	14.8-20.6	18.5	15.7-21.7	26.6	198.4	113



Trunk supports with straps
ITEM # SP010 - SP011
Provides additional trunk support



Pommel
ITEM # PL020 - PL021
Encourages hip abduction



Tray retaining strap set
ITEM # ST056
Secures the tray and keeps the tray in the desired position



Moulded pelvic support
ITEM # SP017 - SP019
Provides snug contoured support



Samba style headrest
ITEM # HT014 - HT015
Provides extra support for the head
Can be shaped by adjusting the sides
Available in two sizes:
Size 1 - 12.7 (w) x 4.7 (h)
Size 2 - 15.2 (w) x 4.7 (h)



Footplate & heel huggies
ITEM # FP017 - FP019
Provides extra foot support

Hardrock chair comes complete with mobile base, 4 point pelvic belt and either a high or low back (8203LB) slotted upholstery as standard.



heathfield

The perfect start

Heathfield is the perfect starter chair for very young children with sizes up to young adults who require minimal levels of postural support when seated. Its unobtrusive design means that it fits well into nursery, school and home environments.

A natural choice

We believe in wood as a natural material for our chairs. Wood is extremely durable and it holds a contemporary beauty that makes it aesthetically pleasing while being environmentally friendly. Our timber comes from sustainable sources and offers thermal properties that make it warm to the touch.

Meeting your child's needs

The Heathfield and its many accessories work in partnership together to provide support and comfort for your child. We recognize that every child is different and the numerous configurations of the Heathfield will ensure that the chair is tailored to your child's needs.



The Heathfield chair comes complete with a buckle pelvic strap, adjustable armrests and a contoured seat.

Heathfield technical data

all dimensions in inches unless otherwise stated

Item #	Size	Seat to floor*	Width	Seat depth	Back height	Max user weight (lbs)	Product weight (lbs)
6901	1	7.1 - 9.1	9.4	7.1 - 9.4	9.8	55.1	11
6902	2	7.9 - 11.8	11.0	8.7 - 12.2	11.0	99.2	13
6903	3	9.4 - 14.6	12.8	10.2 - 13.4	12.6	110.2	18
6904	4	11.4 - 15.4	14.4	12.6 - 15.7	14.4	132.3	35
6905	5	12.8 - 16.7	17.1	14.6 - 18.1	16.5	154.3	44

* Height of chair with contoured seat. The height with a flat seat is 1 inch less

heathfield accessories




Skis

ITEM # SK001 - SK005
Provides greater stability to avoid tippy-ness



Ski & full footrest

ITEM # SK023 - SK027
Provides foot support and offers greater stability



Mobile dolly base

ITEM # DL006 - DL010
Easy to move with swivel multi-directional casters



Reversible footboard

ITEM # FB006 - FB010
Foot support that prevents unwanted chair movement



Adjustable footrest

ITEM # FR005 - FR009
Easy to adjust footrest
The sides of the footrest help to keep feet within the chair



Minimum raise mobile base

ITEM # MB005 - MB007
Facilitates standing transfer
Easy to maneuver
Only available with sizes 3 – 5 chairs



Four point harness

ITEM # HN001 - HN004
Supports the trunk and can provide stability.



Groin harness

ITEM # HN039 - HN041



heathfield accessories




Pelvic cushions – double thickness (pair)¹
ITEM # CS050 - CS053, CS230-
CS234, CS077 - CS080
Decreases seat width
Available in various sizes



Pelvic cushions – single thickness (pair)²
ITEM # CS005 - CS008, CS014 -
CS017, CS023 - CS026
Provides hip and thigh support
Available in various sizes



Pommel⁵
ITEM # PL001 - PL002, PL013 - PL014
Encourages hip abduction



Footplate & straps³
ITEM # FP001 - FP005
Ensures feet are firmly in position



Knee blocks^{3,5}
ITEM # KN001
Maintains knee position
Encourages hip abduction
Only available with sizes 3 - 5 chairs



Adjustable height push handle
ITEM # PH010 - PH014
Adjusts to achieve optimum push position.
To be used in conjunction with a Mobile
dolly base or Minimum raise mobile base



Tray⁴
ITEM # TR001 - TR005
Promotes trunk support
Additionally can be used for activities,
playing and feeding



Quick release tray⁶
ITEM # TR101 - TR104
The tray comes in white and only
available on chairs sizes 1-4

1 Available in sizes 1+1in, 1.5+1in & 2+1in
2 Available in size 1in, 1.5in & 2in
3 Only available in conjunction with a Mobile dolly base,
Reversible footboard or Skis with half & full footrest
4 Heathfield chair sizes 1 & 2 must have a suitable stabilizer
(footrest, mobile dolly base or skis) when attaching trays
5 The Pommel and Knee Blocks cannot be used together
6 Heathfield chair sizes 1 & 2 must have a suitable stabiliser
e.g. footrest/dolly/skis in order to accommodate tray

corner chair



Maximum adjustability
The Corner Chair offers high level of adjustability in the back, side angle, side pad, headrest and pommel. This means that the chair can be positioned to suit your child's individual needs while also enabling the chair to easily grow with your child. Maximum user weight of 187lbs.

Floor level support
Our Corner Chairs are ideal for young children who need to develop sitting balance. During floor sitting some children are inclined to fall back or sideways if they overstretch or lose concentration. Our Corner Chairs are designed to provide support at the back and sides of children to improve their sitting balance. This encourages social interaction and allows them to engage in floor-based activities without the risk of falling over. The pommel at the front is designed to keep legs abducted to promote a stable posture.

Corner Chair technical data
all dimensions in inches

Item #	Size	Height range	Width	Seat depth
5501	1	16.1-20.4	15.7	11.8
5502	2	20.4-25.9	15.7	11.8
5503	3	21.7 - 27.2	17.7	13.7

Support and comfort
Our range of Hip Spica chairs has been specifically designed for young children who are in plaster casts following hip surgery and require supportive seating.

Portable Hip Spica chair
The new design offers high levels of adjustability and creates a lightweight solution. The Portable Hip Spica features a new slimline design and can easily be dismantled so it's simple to transport. It is supplied with a 5-point harness and quick release tray as standard.



Portable Hip Spica chair technical data
all dimensions in inches

Item #	Size	Seat height	Footprint	Tray size	Max user weight lbs	Product weight lbs
5596	1	8.3 - 14.2	22.0x26.0	21.7x17.3	55	22

felix

simple
stackable
school chair



Comfortable learning

Standard classroom chairs can be uncomfortable for a whole day of learning but Felix comes with a seat and back cushion all in one, this gives that little extra padding on the seat and back. The Felix chair provides postural support when needed, when a child is comfortable learning becomes more effective.

Simple features

Felix is stackable, an important feature for any classroom at the end of a working day. The cushions are easily removable and wipeable. They come in six vibrant colors and the chair comes in five sizes. Felix also comes flat packed with easy to follow assembly instructions reducing packaging and making it easier for delivery and handling.

Classroom chair

Felix has been designed as an active classroom arm chair. The lateral side arm rests not only provide support while seated but it also facilitates the sit to standing movement and vice versa. The shape of the seat and back promotes good posture and the design of the chair legs provides stability and prevents tipping.


Felix technical data

all dimensions in inches unless otherwise stated

Size	Seat to floor height	Seat Width	Seat depth	Max user weight (kg)	Product weight (kg)
1	10	11 1/2	8 1/2	165	8
2	12	13	10 1/2	176	10.5
3	14	14 1/2	12 1/2	198	14
4	16	16	14 1/2	242	18
5	17 1/2	16	14 1/2	242	22

Colours	Size 1 Item #	Size 2 Item #	Size 3 Item #	Size 4 Item #	Size 5 Item #
pink	8166	8167	8168	8169	8170
teal	8171	8172	8173	8174	8175
lime	8176	8177	8178	8179	8180
cobalt	8181	8182	8183	8184	8185
red	8101	8102	8103	8104	8105
blue	8161	8162	8163	8164	8165

juni



for dynamic learning

Whether in the classroom, nursery, or at home, good posture when seated creates a positive experience for your child; fine motor skills develop, confidence increases and the ability to engage and learn becomes easier. The Juni chair has a contemporary appeal that helps to promote social inclusion.

Highly adjustable

We understand that your child’s needs change during each stage of growth and development. Seat height and depth as well as armrest height can all be easily adjusted. The backrest can be reclined by 5° increments (from 95° to 105°), while it also features an unique built-in lumbar support.

Posture at its foundation

Juni’s seat and backrest both provide unique contouring. This promotes better posture, gives higher levels of comfort and helps your child feel more secure and stable in the Juni. It is our vibrant chair for children who require minimal levels of postural support when seated.

Dolly Base

ITEM # PH101 - PH104
Facilitates manoeuvrability of your Juni chair
Available in 4 sizes



The Juni chair comes with contoured seat and backrest, built-in lumbar support, pelvic strap, height and depth adjustable seat, height adjustable armrests as well as angle adjustable backrest as standard.

Juni technical data

all dimensions in inches unless otherwise stated

Size	1	2	3	4
Seat to floor height	6.3 - 8.7	7.9 - 11.0	9.4 - 13.4	11.8 - 15.7
Width	9.4	11.0	12.8	14.4
Depth	7.1 - 9.4	8.7 - 11.8	10.2 - 14.2	12.6 - 16.5
Back height	10.2	11.4	12.8	15.7
Back angle (°)	95/100/105	95/100/105	95/100/105	95/100/105
Max user weight (lbs)	66.1	88.2	110.2	132.3
Product weight (lbs)	9	13	17	22

Color	Item #'s			
Blue	16031	16032	16033	16034

juni accessories



Skis
ITEM # SK101 - SK104
Helps to stabilize and prevent chair movements.



Mobile skis & footrest
ITEM # MK101 - MK104
Gives clearance to facilitate standing transfer. Includes a sliding footrest & push handle.



Ski & footboard¹
ITEM # SK111 - SK114
Provides foot support and offers greater stability.



Footrest
ITEM # FR051 - FR054
Height and angle adjustable footrest. Can be attached to the Juni or used separately.



4-point harness
ITEM # HN231 - HN234
Provides additional trunk support.



Quick release tray
ITEM # TR101 - TR104
White quick release tray.



Side support cushions
Provides hip and thigh support
Available in 1 or 2 inch sizes. Sold as a pair.
ITEM # - 1 Inch CS531 & CS532 ITEM # - 2 Inch CS631 & CS632



Pommel
ITEM # PL131 & PL132
Encourages hip abduction. Can be removed with a hand tool.

juni accessories

We've worked in close collaboration with therapists, parents and children to ensure that our products meet and surpass your expectations. Through your valuable feedback the Juni also has companions: highly adjustable Juni Table, Juni Hi & Lo Boosters and Juni Hi Lo Base to overcome height issues at home, nursery & school.

Juni Table

The Juni table offers two key benefits; the legs are adjustable so your child can work or play at the correct height and the table folds, which allows for easy storage. Perfect companion to your favorite Juni chair.



Juni table technical data

all dimensions in inches unless otherwise stated

Size	1	2
Height	12.2 - 17.7	15.7 - 23.6
Width/length	27.6	30.7
Depth	20.9	23.2
Color	Item #	
Blue top	2431	2432

* Castors with brakes are sold as an accessory (Item # JB001); this will increase the height of the Juni Booster by 2.56 inches

Juni Hi Lo Rolling Base

We understand the challenges of the classroom and home environment. Height adjustment is now easy with our Juni Hi Lo Base. The slimline footprint of the base makes this the ideal solution to fit under most classroom desks. So now, you can use your Juni chair together with the Juni Hi Lo Base by simply attaching it to the plate on top of the base; or use the chair independently where height adjustment's are not needed.



Juni Hi Lo Base technical data

all dimensions in inches unless otherwise stated

Size	1		2	
Dimension	Min	Max	Min	Max
Seat to floor height	6.3	22.4	7.9	22.4
Lower leg length to footrest	6.3	8.7	7.9	11.0
Push handle height	19.7	39.4	19.7	39.4
Base length and width	25.5 x 20.5		25.6 x 20.5	
Max user weight (lbs)	66.1		88.2	
Product weight (lbs)	23.1		23.1	
Color	Item #'s			
Natural Footrest	HL141		HL142	

zoomi

Benefits of Zoomi

- Economical seating solution
- Grows with your child
- Highly adjustable
- Tailored to the needs of your child
- Contemporary design
- Easy to clean
- Stable base with slimline footprint
- Can be used at home and in school


Our standard Zoomi

The standard Zoomi is ideally suited to the classroom and home. Both the seat and footrest can be easily adjusted to the appropriate height and depth, allowing the chair to accommodate your child’s needs as they grow.

Zoomi technical data

all dimensions in inches unless otherwise stated

Size		1
Seat to floor height range		4.3 - 24.2
Seat width		17.7
Seat depth		21.7
Back height		10.2
Back angle (°)		95/100/105
Frame height		31.9
Max user weight (lbs)		187.4
Product weight (lbs)		35
Colour	Zoomi standard	Zoomi +
Blue	13031	13032



Standard features

- Frame
- Low backrest
- Flat laminated seat
- Laminated footboard
- Pelvic strap

Perfectly customizable

With multiple options to cater for a variety of users, Zoomi can be uniquely configured to meet the needs of your child. We have put together two simple options:

- Standard Zoomi
- Zoomi+

zoomi+

One size fits all

Zoomi’s innovative and unique seating system has been designed to grow with your child, covering an extensive age range, all the way up to young adults.

Zoomi+

Following collaboration with a number of therapists, we have developed Zoomi+. Complete with low back padded cover, armrest, quick release tray, 5-point harness, flat padded seat, slotted footrest and mobile skis, it offers the ideal solution for children who require a little extra postural support.

Standard features

- Frame
- Low back padded cover
- Armrest
- Quick release tray
- 5-point harness
- Flat padded seat
- Slotted footboard
- Mobile skis




Armrest

ITEM # AR019

Provides arm support

Used for the Quick release tray



Quick release tray

ITEM # TR104

Helps to stabilize trunk

Available in white only

Must be used with armrest

seating • standing • desks & stools • toileting & bathing • changing & showering

colors

The majority of our seating products are typically stocked in blue.

While alternative color choices are available, these options constitute custom products, resulting in extended lead times due to their personalized manufacturing and shipping process (unless otherwise stated in the catalog).



USA Contact Details:

telephone

212 541 2455

email

info@smirthwaite.com

website

www.smirthwaite.com



Smirthwaite USA LLC

405 Lexington Avenue, 26th Floor

New York NY10174



smirthwaite

Together we
transform lives

Smirthwaite reserves the right to modify product specifications at anytime without prior notification.

All measurements are approximations and can vary depending on the wood and materials used.

All details shown are correct at time of print. Smirthwaite® is a registered trade mark of Smirthwaite Limited.