

buddy

Supportive and affordable standing

Smirthwaite's new stander, the Buddy, is an entry-level vertical stander. Designed for younger children (please see technical data for sizing) who require low to moderate levels of support while weight bearing. Its design allows it to be used as a hip stander suitable for targeted therapy for the development of specific areas such as trunk and shoulder girdle stability and head control.

Inclusive and accessible

The Buddy stander is slimline so it can be used in schools, enabling your child to perform numerous activities without feeling isolated from friends. Its color and friendly character also brightens up any room. With the economic cost we are delivering an entry-level stander that is an accessible option to everyone.

Adjustability

The Buddy stander offers numerous height and width adjustments in the chest and pelvic bands as well as additional angle adjustments in the knee supports and sandals. The Buddy size 2 has an additional height adjustment for taller or faster-growing children.

Buddy stander technical data

all dimensions in inches unless otherwise stated

Size	Size 1 Buddy	Size 2 Buddy 2
Item #	21001	21002
Thoracic support (top) to foot board height	18.9 - 24.4	27.6 - 37.0
Pelvic support (top) to foot board height	13.8 - 18.1	19.9 - 29.3
Knee supports (top) to foot board height	7.1 - 13.4	6.3 - 22.0
Foot board to floor height (excluding sandals)	4.7	4.7
Max width (between thoracic supports)	17.7	17.7
Tray size (depth and width)	16.3 x 20.5	16.3 x 20.5
Base frame (length, width and height)	25.6 x 19.6 x 30.9	31.5 x 24.6 x 39.1
Base frame (length, width and height) (high setting)	N/A	31.5 x 24.6 x 43.8
User weight limit (lbs)	88.2	132.3
Product weight (lbs)	30.9	46.3

Buddy and Buddy 2 come complete with tray, chest band, pelvic band, knee supports and sandals. Straps and tray only available in lime.




star

The Star, a vertical and prone stander

The physiological benefits of standing are enormous – not to mention the psychological benefits of being able to tackle activities. The Stander aims to make standing enjoyable and practical for your child, while making it simple to achieve.

The Stander is a vertical and prone entry-level stander that enables children from five years old to young adults to gradually develop and maintain their ability to weight bear by promoting excellent symmetry. It is suitable for children who require lower to higher levels of postural support and control. The Stander is ideal for both school and home environments.

The Stander offers high levels of adjustment. Anterior, pelvic and knee supports are height,



Star technical data

all dimensions in mm unless otherwise stated

Size	Size 1 Star
Code	21003
Thoracic support (top) to foot board height	27.2-39.3
Pelvic support (top) to foot board height	18.8-31.8
Knee supports (top) to foot board height	8.7-13.5
Foot board to floor height (excluding sandals)	7.5
Max width (between thoracic supports)	18.1
Tray size (depth and width)	15.3 x 17.5
Base frame (length, width and height)	1006 x 660 x 1200
Base frame (length, width and height) (high setting)	39.3 x 25.9 x 36.2/45.6
User weight limit (lbs)	60
Product weight (lbs)	28

Tilt angle on Star: 40° - upright.

