



Toileting & Bathing



standing • seating • desks • stools
toileting • bathing
changing • showering



supporting children and
their families for life

bath corner chair



The Bath Corner Chair comes complete with velcro pelvic strap, adjustable pommel, telescopic back rail and telescopic handrail as standard.

More support

Our Bath Corner Chair is designed for children who require a little extra help with sitting balance when bathing. This is achieved by its angled back and hand rail, both of which provide additional security.

Adjustable

Like the Bath Chair, our Bath Corner Chair provides numerous levels of adjustment. The pommel and pelvic strap can be adjusted to suit each individual child's needs.


Compact

Lightweight and easy to clean, the Bath Corner Chair has a perfectly compact design. This allows it to sit low in the bath enabling your child to fully extend their legs in complete comfort. Its size also enables it to fit a wide range of baths.

Bath Corner Chair technical data

all dimensions in inches unless otherwise stated

Size	1	2
Item #	7511	7512
Seat height	3.0	3.0
Seat width	13.4	13.4
Seat depth	11.0	13.0
Back height	13.8	18.1
4 Pommel settings (distance to back of pommel)	4.3 / 5.1 / 6.3 / 7.5	5.5 / 6.7 / 7.9 / 9.1
Max user weight (lbs)	66.1	77.2
Product weight (lbs)	6.6	8.8



combi



Progressive toileting

Our Combi toileting chair is a great starter chair. Its design enables it to be used as a supportive chair while at the same time providing numerous features to assist your child with the everyday toileting needs.

Providing solutions

The Combi toileting chair's design also incorporates a potty for younger children. This can easily be removed to provide a commode seat for children who have progressed their toilet training.

Adjustability

Good levels of adjustability are offered with the Combi toileting chair so that a wide range of needs can be met. Seat height and depth can be adjusted to suit the needs of your child.

The Combi Chair comes complete with potty, plastic commode seat, upholstered seat and velcro pelvic strap as standard.

Combi technical data

all dimensions in inches unless otherwise stated

Size	1	2	3	4
Item #	7551	7552	7553	7554
Seat to floor height*	8.5 - 11.4	10.8 - 14.8	11.8 - 15.7	13.4 - 17.7
Seat width	11.4	12.8	14.4	17.1
Seat depth	8.1 - 11.4	9.6 - 14.0	12.2 - 15.2	13.8 - 17.3
Back height	12.8	14.6	15.7	19.7
Max user weight (lbs)	77.2	99.2	121.3	143.3
Product weight (lbs)	24.3	28.7	30.9	33.1

*Seat height to floor is to the top of upholstered seat (deduct 1in for the height of the Potty and deduct 1.6in for the height of the Commode seat)

combi accessories



Skis & full footrest¹
 ITEM # SK037 - SK040
 Provides foot support and offers greater stability



Reversible footboard¹
 ITEM # FB023 - FB026
 Provides foot support
 Prevents unwanted chair movement



Independent adjustable footrest¹ ITEM # 4187
 Portable, height and angle adjustable
 Solid support for the feet



Skis
 ITEM # SK051 - SK054
 Provides greater stability




Footrest & straps¹
 ITEM # FP002 - FP005
 Ensures feet are firmly in position



Splash guard
 ITEM # SG001 - SG002
 Assists in leg abduction and keeps toileting contained
 Requires Toggle handrail



Pelvic cushions single or double thickness (pair)²
 Provides hip and thigh support
 Available in various sizes



Toggle handrail
 ITEM # HR014 - HR017
 Allows child to hold and maintain balance and to sit with minimal support



Tray
 ITEM # TR003 - TR006
 Provides trunk support
 Additionally can be used for activities such as playing

¹ Only available in conjunction with a Mobile dolly base, Reversible footboard or Skis with half & full footrest
² Available as double cushions (pair) size 1+1in, 1.5+1in & 2+1in and singular cushions (pair) size 1in, 1.5in & 2in - please ask for ITEM #'s

seating • standing • desks & stools • toileting & bathing • changing & showering

a helping hand around
the World



always here to help

USA Contact Details:

telephone

212 541 2455

email

info@smirthwaite.com

website

www.smirthwaite.com

Smirthwaite USA LLC

405 Lexington Avenue, 26th Floor

New York NY10174



Smirthwaite reserves the right to modify product specifications at anytime without prior notification.
All measurements are approximations and can vary depending on the wood and materials used.
All details shown are correct at time of print. Smirthwaite® is a registered trade mark of Smirthwaite Limited.