

**smirthwaite**

Together we  
transform lives



# Essential changing bench

A cost-effective,  
functional pediatric  
changing bench



[smirthwaite.com](http://smirthwaite.com)

# Essential Changing Bench Overview

smirthwaite

Introducing the Essential fixed-height changing bench to support dealers seeking a competitively priced solution. This model provides a high-quality, affordable alternative within the pediatric market.



## Product features

- Superb low transfer height for users
- Excellent height range from 18.1 in to 38.98 in
- Safe working load and lifting capacity of 400.92 lbs (200kg)
- Choice of 2 sizes: 59.8 in or 72.4 in
- 3.94 in diameter individual locking castors
- Side support rails fitted as standard on both models
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard

## Product dimensions

### 59.8 in / 152 cm Model

Length:	59.8 in (152 cm)
Width:	25.6 in (65 cm)
Height:	18.1 in (46 cm)
Safe Working Load:	440 lbs (200 kg)
Safety:	Includes padded safety arms

### 72.4 in / 184 cm Model

Length:	72.4 in (184 cm)
Width:	25.6 in (65 cm)
Height:	18.1 in (46 cm)
Safe Working Load:	440 lbs (200 kg)
Safety:	Includes padded safety arms

## Colour

Standard vinyl in Canard.



# Why choose Essential?



## High-quality & Durable

Robust construction designed for reliability in care environments.



## Easy to Clean

Simple to wipe, long-lasting vinyl.



## Ergonomic & Safe

Adjustable height supports caregiver ergonomics. Safety rails and high load capacity provide user security.



## Mobile & Flexible

Ideal for schools, clinics, changing facilities, and other care environments.



## Cost Competitive

Cost-effective option without compromising essential function.



## Long Life Span

12 years with proper maintenance and normal use.



# Our advantage

Smirthwaite gains several key advantages that reinforce our value proposition



## Scale & Stability

We are a well-established, trusted paediatric company, ensuring consistent supply, strong infrastructure, and guaranteed aftercare.

## Shared Clinical Knowledge

We benefit from cross-brand expertise in seating, moving & handling, and bathing, helping us deliver clinically integrated solutions.



## Investment in Quality & Innovation

We invest in better materials, smarter modular design, and more robust support systems, resulting in longer-lasting equipment and lower running costs.

## Unified Ethical Vision

Our company-wide commitment to sustainability ensures our products are designed for recycling, reuse, and adaptability, aligning with local government goals.



Together we  
transform lives

## Contact

212 541 2455

info@smirthwaite.com

smirthwaite.com



**smirthwaite**

Together we  
transform lives