

Essentials Changing Bench

The Essentials Changing Bench provide a variable height support surface to reduce back strain and posture problems for carers and support staff during routine changing procedures in environments such as special schools, specialist adult centres and paediatric units. These tables are available in two sizes and supplied with side support rails as standard. They are fitted with large wheels with brakes for ease of transport and positioning.



Height range 18.1 in to 38.98 in




63

59.8



63

72.4



400.92
lbs



Colour: Canard



Optional Side Rail Cushions can be easily removed from the Side Rails for cleaning

Standard features

- Safe working load and lifting capacity of 400.92 lbs
- Excellent height range from 18.1 in to 38.98 in
- Superb low transfer height for patients
- Choice of 2 lengths, 59.8 in or 72.4 in
- Table width - 25.6 inches
- 3.94 in diameter individually locking castors
- Side support rails fitted as standard on all models
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on all models

Installation

The Table is designed as a "plug-and-play" system to ensure rapid deployment and ease of use. The unit is delivered fully assembled with all primary mechanical and electrical components pre-commissioned.

Maintenance Note: As the side cushions are fitted using hand wheels, they can be easily removed without tools for deep cleaning or disinfection to maintain clinical hygiene standards.