



Therapy & Stools



seating • standing • toileting
bathing • changing • showering
desks • tables • therapy

smirthwaite

Together we
transform lives

welcome

Since 1986, Smirthwaite has developed and retained its reputation as one of the leading designers & manufacturers of pediatric equipment. In that time, we have dedicated ourselves to finding and providing effective solutions that offer significant benefits to babies, children and young adults in their daily lives. Working closely with children, parents, carers and therapists, we have developed a comprehensive range of seating, standing, tables, desks, toileting, bathing, changing, showering and therapy products.



Innovation in our design

Special needs equipment is often overly complex, difficult to use and expensive to maintain. At Smirthwaite, we believe in well-made, easy to use products that stand the test of time. Through close consultation with leading therapists, parents and end users, we ensure our designs offer numerous clinical benefits whilst maintaining a customer-centric and innovative approach.

Easy to clean

All of our products are easy to clean. For instance, our wooden chairs have anti-microbial properties that ensure surfaces are kept hygienic. This helps to reduce any risk of cross infection and is especially useful in multi-user environments such as nurseries and schools.

All our upholstery is carefully chosen due to its clever material make up that contains antibacterial agents that help against the spread of MRSA, E coli and other bacteria.

Warranty & servicing

All of our products come with a two year warranty as standard. An inspection, maintenance, service and repair facility for all Smirthwaite products is available to you, for complete peace of mind. This can be scheduled at a time that is convenient for you.

Call our Customer Services team today on 212 541 2455 or e-mail info@smirthwaite.com to book your service, inspection or repair.

Therapy & stools

Therapy plays a central role in helping children to keep active, manage symptoms and recover from illness or injury.

We can provide therapy equipment suitable for many different types of therapy, helping children to unlock their potential and maximise their development. Our equipment provides fun, comfort and achievement.

contents

therapy & stools

footrests, steps & bars	5
therapy bench & bolster	6
wheeled stools	7
accessories	8





footrests, steps & bars



Interlocking steps
CODE 4188

Interlocking steps

Our 3 Interlocking Steps come complete with a single handrail that can be locked into position. All 3 steps have the same width 24.61 inches and depth 14.57 inches. The lower step is at a height of 4.92 inches, the middle step at 9.84 inches and the platform at 14.76 inches. Max user weight of 220.46 pounds and product weight of 44.09 pounds.

An extra handrail (CODE 4189) for the Interlocking steps is available to buy separately.

improving your child's balance

Parallel bars

CODE 4292 - 4293

Our Parallel Bars are height and width adjustable and aim to improve your child's balance. They are extremely easy to use and adjust.



Parallel Bars technical data

all dimensions in mm unless otherwise stated

Size	1	2
Code	4292	4293
Length (inches)	78.74	104.72
Height min (inches)	14.17	19.69
Height max (inches)	21.26	32.28
Width min (inches)	14.96	19.69
Width max (inches)	23.62	28.74
Max user weight (lbs)	176.37	220.46
Product weight (lbs)	220.46	55.12

Interlocking footrest

CODE 4185-4186

Our Interlocking Footrest is portable and can be used with a number of seating or therapy products. There are three interlocking footrests and they can be used together or separately. Available in 2 sizes.



Interlocking Footrest technical data

all dimensions in mm unless otherwise stated

Size	1	2
Code (inches)	4185	4186
Width (inches)	15.35	16.54
Depth (inches)	13.39	11.81
Height min (inches)	1.97	2.95
Height max (inches)	5.51	9.06
Max user weight (lbs)	176.37	176.37
Product weight (lbs)	11.02	11.02



therapy bench

“For our children, we look for products that will have great clinical benefits for the management of their posture.”

Anna, Teacher, Midlands, UK

Our popular height and angle adjustable Therapy Bench together with accessories¹ is designed to assist in a wide range of children’s therapy programmes both at home and school.

Therapy Bench technical data

all dimensions in inches unless otherwise stated

Size	1	2	3	4	5	6	7	8
Code	5541	5542	5543	5544	5545	5546	5540	5547
Length	24.02	24.02	32.28	32.28	32.28	32.28	26.77	31.50
Width	9.84	11.81	11.81	12.80	14.96	14.96	11.02	14.17
Min height	7.68	7.68	11.02	11.02	15.94	15.94	6.69	10.63
Max height	11.61	17.91	27.76	27.76	27.76	27.76	10.63	15.75
Max user weight (lbs)	220.46	220.46	220.46	220.46	220.46	220.46	220.46	220.46
Product weight (lbs)	15.43	17.64	24.25	26.46	35.27	22.05	22.05	24.25



Colours



therapy bolster



Therapy bolster

CODE 5861 - 5866



Low end (single)

CODE LG014
(W 500mm x H 450mm)



High end (single)

CODE LG013
(W 500mm x H 700mm)

Colours



Our Therapy Bolster is extremely versatile and can be used just as a roll or built up to have a combination of High ends or Low ends.

Therapy Bolster technical data

all dimensions in inches unless otherwise stated

Size	1	2	3	4	5	6
Code	5861	5862	5863	5864	5865	5866
Length	23.62	47.24	23.62	47.24	23.62	47.24
Height min	11.81	11.81	12.80	12.80	13.78	13.78
Height max	18.70	18.70	19.96	19.96	20.67	20.67
Diameter	7.87	7.87	9.84	9.84	11.81	11.81
Max user weight (kg)	220.46	220.46	220.46	220.46	220.46	220.46
Product weight (kg)	13.23	24.45	17.64	28.66	22.05	330.07

¹Therapy Bench accessories, as shown in the image on the page opposite: Adjustable backrest with pelvic guides CODE BK091 and Adjustable femoral kneeblocks (pair) CODE KNO02

wheely stools



adjustable height wheely stool

Colours

- Red upholstery & yellow frame
- Green upholstery & blue frame
- Pink upholstery & blue frame



round wheely stool

Colours



Space saving for therapists

Our space saving, Round Wheely Stool is a popular choice for therapists, teachers and parents when working at a low level with children. The stool comes in our natural wood finish and is available in one size only.

Therapist’s favourite

The Adjustable Height Stool is our most popular wheely stool. It has been designed specifically for working at lower levels with your child, student or client. The stool is designed for sitting astride the narrower side. Available in 3 colour options with a height range of 13.39 - 18.11 inches. One size only.

Adjustable Height Wheely Stool technical data

all dimensions in inches

Code	Size	Colour	Width	Height min	Height max	Max user weight (lbs)	Product weight (lbs)
4181	1	Red & yellow	17.7	13.4	18.1	220	15.43
4182	1	Green & blue	17.7	13.4	18.1	220	15.43
4184	1	Pink & blue	17.7	13.4	18.1	220	15.43
4192	1	Red & yellow	17.7	13.4	18.1	220	17.64
4193	1	Green & blue	17.7	13.4	18.1	220	17.64
4194	1	Pink & blue	17.7	13.4	18.1	220	17.64

Round Wheely Stool technical data

all dimensions in inches

Code	Size	Colour	Height	Diameter	Base	Max user weight (lbs)	Product weight (lbs)
4178	1	Pacific blue	15.75	14.17	18.90 x 18.90	220.26	19.84
4179	1	Deep pink	15.75	14.17	18.90 x 18.90	220.26	19.84



seating standing toileting bathing
changing showering desks tables therapy


accessories



Suction hand grabs

CODE 4201-4203

Used to stabilise the hand to increase function in the opposing hand



Dycem non-slip mats

CODE 4703 red or blue (250 x 180mm)

CODE 4706 red or blue (350 x 250mm)

CODE 4709 red or blue (450 x 380mm)

CODE 4700 red or blue (round 190mm)



US Contact Details:

telephone

+ 212 541 2455

email

info@smirthwaite.com

website

www.smirthwaite.com

Smirthwaite Ltd

16 Wentworth Road, Heathfield,
Newton Abbot, Devon TQ12 6TL

smirthwaite

Together we
transform lives



**MADE IN
BRITAIN**



Smirthwaite reserves the right to modify product specifications at anytime without prior notification. All measurements are approximations and can vary depending on the wood and materials used. All details shown are correct at time of print. Smirthwaite® is a registered trade mark of Smirthwaite Limited.